Roche CARDIAC IQC

cobas®
SYSTEM
cobas h 232

REF

04880668190

1 Roche CARDIAC IQC high test strip
1 Roche CARDIAC IQC low test strip

English

Intended use

The Roche CARDIAC IQC is an *in vitro* diagnostic device and serves as a performance check for the optical system of the cobas h 232 instrument.

The Roche CARDIAC IQC is intended for near-patient testing.

Not for self-testing.

Summary

The Roche CARDIAC IQC test consists of two re-usable control strips to ensure correct functioning of the instrument.

The Roche CARDIAC IQC low test serves to check the lower measurement range.

The Roche CARDIAC IQC high test serves to check the higher measurement range.

Performance checks of the optical system should ideally be carried out on a daily basis. The cobas h 232 instrument should be checked using the Roche CARDIAC IQC test according to the applicable local or national regulations.

Reagents

One test contains:
Gold-labeled anti-HSA antibodies
IQC low: 85 ng
IQC high: 126 ng

Precautions and warnings

For *in vitro* diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

- Do not touch or wipe the signal line area of the test strip.
- Do not apply any sample material to the test strip.
- Do not expose the test strip to sunlight.
- Protect the test strip from dust and sunlight.
- Close the container immediately after removing the test strip to protect the second test strip from dust and moisture.

▪ Remove the test strip (low or high) from the cobas h 232 instrument directly after the measurement is performed and place it quickly into the container to protect it from dust and moisture.

▪ For information on the optimum ambient temperature for the cobas h 232 instrument, consult the operator's manual for the instrument.

Lot code

Every kit contains a lot-specific code chip. The instrument display prompts the user to insert the chip. To ensure that the code chip and test strip lot match, compare the lot number in the display with the number on the code chip. The code chip provides the cobas h 232 instrument with all the required lot-specific information. An error message is displayed if the wrong code chip is inserted for a test strip lot.

Test procedure

1. Bring the Roche CARDIAC IQC test strips to room temperature before starting the measurement.
2. Switch on the cobas h 232 instrument.
3. Remove a test strip, like Roche CARDIAC IQC low, from the container.
4. Close the container immediately with the desiccant stopper.
5. Choose "QC Test" from the main menu.
6. When instructed in the display, insert the test strip with a smooth movement while holding it in an horizontal position. Be careful not to scratch the barcode on the underside of the strip. This can also happen if you push the strip over the small plastic positioning bar during insertion.
7. Measurement starts as soon as the test reaches the measuring position.
8. The result is displayed after a few seconds.
9. Place the measured test strip quickly into the container and repeat the measurement with the second Roche CARDIAC IQC test strip from the container (start from step 3).

For more details on the procedure, refer to the operator's manual.

Storage and stability

Unopened, at 2-30 °C up to the stated expiration date

After opening 6 months

Write the opening date on the container label.

Materials provided

- REF 04880668190, Roche CARDIAC IQC
- 1 Roche CARDIAC IQC high test strip
- 1 Roche CARDIAC IQC low test strip
- 1 code chip

Materials required (but not provided)

- REF 04901126xxx, cobas h 232 instrument
- REF 04901142xxx, cobas h 232 instrument with scanner

Display of results

On completion of measurement, the cobas h 232 instrument shows the result in the display. The time required for the Roche CARDIAC IQC test to display a value is approximately 20 seconds. The cobas h 232 instrument determines the reflectance value of the signal and displays "Pass" or "Fail".

"Pass" means that the optical system of the cobas h 232 instrument is functioning properly in the tested measurement range. If "Fail" appears, the value is outside the confidence interval.

Values outside the confidence interval:

- Repeat the instrument check using a new Roche CARDIAC IQC test strip of the same type (low or high) from a new, unopened container.
- If the results obtained with the new test strip from the new container fall within the confidence interval, the first test strip is termed defective and must be discarded.
- If the problem persists, contact your local Roche Diagnostics Technical Support team. For further information, please refer to the appropriate Operators Manual for the analyzer concerned, and the Method Sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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If the results obtained with the new test strip from the new container fall within the confidence interval, the first test strip is termed defective and must be discarded.

If the problem persists, contact your local Roche Diagnostics Technical Support team. For further information, please refer to the appropriate Operators Manual for the analyzer concerned, and the Method Sheets of all necessary components.

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For completion of measurement, the cobas h 232 instrument shows the result in the display. The time required for the Roche CARDIAC IQC test to display a value is approximately 20 seconds. The cobas h 232 instrument determines the reflectance value of the signal and displays "Pass" or "Fail".

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