




SMARTsat®

Mehrweg Pulsoximetriesensoren

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en Instructions For Use

SoftCap®

Type: **SC7500**

REF 6020132004

Type: **SCM7500**

REF 6020132010

Type: **SCP7500**

REF 6020132300

SoftWrap®

Type: **W7500**


REF 6020132006

SoftFlap®


Type: **SF7500**

REF 6020132002


(A):


SC7500, SCM7500 (>20kg)


SCP7500 (10 - 20kg)


SF7500 (>20kg)

W7500 (>20kg)


W7500 (10-20kg)


W7500 (<3kg)

(B)


SC7500, SCM7500, SF7500


SCP7500, W7500

(C)


SC7500, SCM7500, SCP7500

(D)

W7500

(E)

SC7500, SCM7500, SCP7500

(F)

SC7500, SCM7500, SCP7500

SMARTsat®

Mehrweg-Pulsoximetriesensoren

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Reusable Pulse Oximetry Sensors

Directions for Use

Note: This Instruction for Use is intended for the medical professional. When prescribing to lay users, ensure that they have been trained regarding all covered warnings, cautions and instructions for use and have a home use guide. For home use, a lay operator instruction for use with is available.

Indications for Use

The SMARTsat® Reusable Pulse Oximetry Sensors SoftCap® SoftFlap® and SoftWrap® are intended for use with monitoring systems that use SMARTsat® technology.

The SMARTsat® Reusable Pulse Oximetry Sensors are intended for use in non-invasive continuous or spot-check monitoring of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate.

The use environment may include: hospitals, hospital-type facilities, intra-hospital transport, mobile emergency medical applications including ground transport and home use.

Transport environments include intra-hospital transport and ground emergency transport (road ambulances). The use environment is dependent upon the use environments specified for the monitoring system utilized in conjunction with the sensor.

The SMARTsat® Reusable Pulse Oximetry Sensors are for prescription use only.

Disposal

Do not dispose of with domestic waste. Follow local regulations when disposing or recycling this device and its component parts. A list of component materials is available by contacting bluepoint MEDICAL Technical Service Department (See cover back page).

Instructions for Use

Intended for use with SMARTsat® technology compatible monitoring systems.

Sensor application**1. Apply SC7500, SCM7500, SCP7500, SF7500:**

Insert the patient's finger into the sensor. The preferred finger to use is the index finger (A). Alternative sites may be used depending on finger sizes (B).

Apply W7500:

Fix a Wrap-Tape to the sensor (D) and apply the sensor to the selected sensor site. The preferred sensor sites are the index finger for adults and pediatrics, the hand for infants and on the foot below the toes for neonates.

2. The finger is correctly inserted when

SC7500, SCM7500, SCP7500, SF7500:

- a. The tip of the finger touches the sensor's inside end. (C)
- b. The sensor cable extends along the top of the patient's hand. (A)

NOTE: Ensure the sensor is not incorrectly installed backwards (E) or with the finger extending (F)

W7500:

- a. The emitter and detector are positioned directly opposite each other
- b. The sensor is applied firmly, but not too tight

3. Then switch on the oximeter and verify proper operation. (ensure oximeter is providing a reading).
4. Visually monitor the sensor site at least every 6 hours to ensure integrity of the skin and that no pressure injury is developing. Reposition to an alternate location if skin integrity changes.

Sensor removal:

1. Open the sensor by pressing on sides of the sensor (SC7500, SCM7500, SCP7500) or at the front part of the housing (SF7500). Remove the SoftWrap® by opening the velcro fastener (W7500). Remove the sensor from the finger and follow the cleaning instructions, if necessary.
2. Disconnect the sensor from the monitor by unplugging the sensors cable.
3. Store the sensor away for its next use.

Sensor expected life is 2 year.

WARNINGS

WARNING: Pulse oximeter probes are designed for use with specific monitoring systems. Incompatible components can result in degraded performance and electromagnetic incompatibility, or patient injury. The operator is responsible for checking compatibility prior to use.

WARNING: A functional tester cannot be used to access the accuracy of pulse oximeter probes.

WARNING: Check the patient skin condition at least every 6 hours and reposition to an alternate location if skin integrity changes.

WARNING: Misapplication of a pulse oximeter probe resulting in excessive pressure (e.g. sensor too small or applied too tight) for prolonged periods can cause pressure injury.

WARNING: Avoid application of the sensors to area with poor skin integrity.

WARNING: Do not sterilize by irradiation, steam, or ethylene oxide – refer to cleaning and disinfection instructions. Use of agents other than specified may damage the sensor.

WARNING: Do not use the sensor if it is damaged. Use of a damaged sensor could cause patient injury or equipment failure.

WARNING: Excessive patient motion, excessive ambient light, electromagnetic interference (e.g. stacking and location to other medical equipment), dysfunctional hemoglobin, intravascular dyes, finger nail polish and long or artificial finger nails may affect the sensor performance and the accuracy of the measurement.

WARNING: Route all cables carefully to reduce the possibility of patient entanglement or strangulation.

WARNING: Do not alter or modify the sensor. Alterations or modifications may affect performance or accuracy.

WARNING: Do not use the SMARTsat® sensor or other oximetry sensors during MRI scanning. Conducted current may cause burns.

CAUTIONS

CAUTION: Avoid possible interference from sources of electromagnetic interference such as, but not limited to: Cellular phones, radio transmitters, motors, telephones, lamps, electrosurgical units, defibrillators, and other devices. Interference may cause inaccurate measurements. If you are unsure if your device is working properly, contact your clinician.

CAUTION: Sensor may work with different monitors. Pay attention for specific monitor IFU for essential performance, accessories, compliance information (EMC) and electrosurgical applications.

Cleaning

It is necessary to clean and disinfect the sensor before attaching it to a new patient. Routine cleaning or low-level disinfection is recommended when using the sensor on the same patient.

Cleaning SC7500, SCM7500, SCP7500, W7500:

To clean the sensor - 70% Isopropyl alcohol:

1. Disconnect the sensor from the monitor before cleaning or disinfecting.
2. Rinse the sensor and cable with distilled water for at least two minutes.
CAUTION: Do not immerse the connector at the sensor cable in liquid.
3. Wipe the cable with a cotton pad soaked in 70% Isopropyl alcohol.
4. Immerse the sensor housing in 70% Isopropyl alcohol for 5 minutes.
5. Allow ten minutes for the sensor, cable, and connector to dry after rinsing.
6. Visually inspect the sensor and repeat the cleaning as needed until a visually clean condition has been reached.

To disinfect (low-level) - 1:10 bleach dilution:

1. Disconnect the sensor from the monitor before cleaning or disinfecting.
2. Rinse the sensor and cable with distilled water for at least two minutes.
CAUTION: Do not immerse the connector at the sensor cable in liquid.
3. Prepare 1:10 bleach dilution in distilled water. Wipe the cable with a cotton pad soaked with the solution of diluted bleach.
4. Immerse the sensor housing in the solution of diluted bleach for five minutes.
5. Rinse the sensor and cable by immersing them in distilled water for five minutes.
6. Allow ten minutes for the sensor, cable, and connector to dry after rinsing.
7. Visually inspect the sensor and repeat the cleaning as needed until a visually clean condition has been reached.

To disinfect (high-level) - CIDEX OPA²:

1. Rinse sensor to remove surface dirt – the sensor housing may be fully immersed in liquid.
CAUTION: Do not immerse the connector at the sensor cable in liquid.
2. Clean the sensor and patient contact surfaces. Prolystica¹ is recommended as a cleaning agent. Follow the manufacturer's instructions for use. Immerse the sensor housing in the solution and wipe the inner and outer surfaces with a soft brush or cloth to remove any visible soil. Rinse sensor with water following cleaning and wipe dry prior to disinfection.
Visually inspect the sensor and repeat the cleaning as needed until a visually clean condition has been reached.
3. CIDEX OPA² is recommended for high-level disinfection. Follow the manufacturer's instructions for use. Immerse the sensor housing in the solution for 12 minutes and remove.
4. Thoroughly rinse the sensor with water following disinfection (5 minutes minimum distilled water rinse/soak).
5. Manually dry the sensor and verify the sensor, cable, and connector are dry prior to use.

1) Prolystica is a registered trademark of Steris Corporation

2) CIDEX OPA is a registered trademark of Johnson and Johnson Corporation

Cleaning SF7500:

To clean the sensor - 70% Isopropyl alcohol:

1. Disconnect the sensor from the monitor before cleaning or disinfecting.
2. Wipe the cable and sensor head at the inside and outside with a cotton pad soaked in 70% Isopropyl alcohol.
3. Allow ten minutes for the sensor, cable, and connector to dry after cleaning.
4. Visually inspect the sensor and repeat the cleaning as needed until a visually clean condition has been reached.

Specifications	LED wavelength	output power³
	Red: 660nm Peak	< 5 mW
	Infrared: 905nm Centroid	< 5 mW
Accuracy	SpO₂ Accuracy:⁴	Pulse Rate Accuracy:⁵
	60 – 100%: Arms ≤ 2	20 – 300 BPM: ±2 BPM (refer to monitor's operating manual for more information)
	60 – 70%: Arms ≤ 3	
	70 – 80%: Arms ≤ 2	
	80 – 90%: Arms ≤ 2	
	90 – 100%: Arms ≤ 2	
	<60%: unspecified	
Water spray protection	SC7500, SCM7500, SCP7500, W7500	
	SF7500	IP22 Protection against diagonal water drips (up to a 15° angle)

Note: Because pulse oximeter measurements are statistically distributed, only about two-thirds of measurements can be expected to fall within ± Arms of the measured values by a co-oximeter. To verify the function of pulse oximeter probes a functional tester like Index 2 or equivalent can be used.

3) **Note:** This information can be especially useful to clinicians

4) Validated by clinical tests, where the measured values of the sensors were compared with those of the arterial co-oximetry in healthy adult subjects over the specified functional oxygen saturation range. These values were validated on the following monitor: OxyTrue®A. A functional tester may not be used to validate SPO₂ accuracy.

5) Pulse rate accuracy was verified by simulated bench tests to ensure entire range was verified. These values were verified on the following monitor: OxyTrue®A

Biocompatibility

The materials that come into contact with the patient have undergone extensive biocompatibility testing. Further information is available on request.

Environmental Specification

Operational temperature range:	<i>SC7500, SCM7500, SCP7500, W7500:</i> +0°C to +40°C (32°F to 104°F) <i>SF7500:</i> +5°C to +40°C (41°F to 104°F)
Transient operating conditions:	Sensor will operate for 20 minutes at -20°C and +50°C temperature conditions.
Storage, and transport temperature range:	<i>SC7500, SCM7500, SCP7500, W7500:</i> -40°C to +70°C (-40°F to 158°F) <i>SF7500:</i> -20°C to +70°C (-4°F to 158°F)

Relative humidity range:	15 - 95% non-condensing
Atmospheric pressure range:	620 hPa to 1060 hPa
Stabilization Time (from storage to operational conditions):	Up to 20 minutes

Note: These environmental specifications are specific to the sensor. For information on the environmental specification of the compatible monitoring systems, please see the Operator's Manual for that system.

Warranty

To obtain information about a warranty, if any, for this product, contact bluepoint MEDICAL's Technical Services Department, or your local bluepoint MEDICAL representative (See cover back page).

Do not use if package is opened or damaged

Manufacturer

Catalog number

Lot number

Do not throw away

Date of manufacture

Serial number

Not made with natural rubber latex

Prescription use only

Not made with DEHP


Consult instructions for use


Follow instructions for use. Symbol appears blue on device.


MR Unsafe


620 hPa
1060 hPa
Atmospheric pressure limitation

15%
95%
Humidity limitation

32°F
0°C
104°F
40°C
Temperature limit

-40°F
-40°C
158°F
70°C
Store at room temperature

Keep away from sunlight

Keep dry

Waste Electrical and Electronic Equipment

0482

Technical Service Department:

For technical information and assistance, contact your local bluepoint MEDICAL representative

V3.2
© 2017 bluepoint medical GmbH & Co. KG.

Made in Germany.

 bluepoint medical GmbH & Co. KG
An der Trave 15
D-23923 Selmsdorf, Germany
www.bluepoint-medical.com