
WESTMED

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Please read this instruction manual for use carefully and completely before use and observe the care instructions.

### **What is blood pressure?**

Blood pressure is the pressure exerted by the blood against the arterial walls. Blood pressure in the arteries is subject to constant fluctuations during a cardiac cycle. The highest pressure in this cycle is called systolic blood pressure, while the lowest pressure is called diastolic blood pressure. You can use these two blood pressure readings to determine your overall blood pressure. Blood pressure can be affected by many different factors, such as physical exertion, anxiety, or the time of day. Blood pressure is usually low in the morning and rises in the afternoon through the evening. It is also lower in the summer than in the winter.

### **What is normal blood pressure?**

Blood pressure varies from person to person and depends on a variety of factors such as age, weight, physical condition, and gender. The classic normal reading for 18- to 45-year-olds is 120/80. Please note: only your doctor is qualified to determine whether the readings you have achieved are normal for you.


## How to measure blood pressure

### 1. Relax in a comfortable position.

- It is very important that the person whose blood pressure is being measured is relaxed and sitting comfortably.
- Measure your arm at heart level with the arm spread out, but not too tense or stiff. Place the arm on a flat surface, such as a table, with the left palm facing up. Keep the arm still during the measurement.

### 2. Positioning of the Velcro cuff and the chest piece

- Place the cuff 2–3 cm above the bend of the elbow on the exposed left upper arm; close the cuff with the Velcro. The cuff should be comfortable but not too tight. One or two fingers should fit between the cuff and the arm. Never place the cuff over clothing.
- For cuffs with a metal D-shaped ring, pull the free end of the cuff through the D-ring and close the cuff with the Velcro.
- For devices for private use, the membrane from the chest piece of the integrated stethoscope should be placed over the brachial artery.

### 3. Inflating the cuff

- Close the air valve on the ball by turning the screw clockwise.
- Press the inflation ball at even intervals until the needle or measuring point is 30 mmHg above your normal systolic pressure value. If you do not know the value, pump up only to 200 mmHg.

### 4. Systolic blood pressure - upper reading.

- Slowly open the air valve by turning the screw counterclockwise and hold the stethoscope chest piece above the brachial artery. A proper air release rate is essential for accurate measurement, so practice and master the recommended release rate of 2-3 mmHg per second, or a drop of one to two marks on the pressure indicator per heartbeat. You should not leave the cuff inflated for longer than necessary. When the cuff begins to deflate, you must listen carefully with the stethoscope. Note the reading on the display as soon as you hear a faint, rhythmic tapping or knocking sound. This is the reading for systolic high blood pressure. Listen carefully and become familiar with the pulse sound. When you get used to it, review the procedure again with your doctor.

## 5. Diastolic blood pressure - lower reading

- Deflate the pressure at a steady rate. When you reach the diastolic low blood pressure reading, the tapping sound will stop.
- Drain the cuff valve completely. Remove the cuff from your arm and the stethoscope from your ears.

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info@fleischhacker.biz

HONSUN (NANTONG) Co., Ltd

No. 8 Tongxing Road Economic & Technological  
Development Area, 226009 Nantong City,

Jiangsu Province, P.R. China

Shanghai International Holding Corp. GmbH (Europe)

Eiffestraße 80, 20537 Hamburg, Germany

SRN: DE-AR-000000001

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Fleischhacker GmbH & Co. KG