

VBM Medizintechnik Entwicklung <i>R&D</i>	Technisches Datenblatt <i>Technical Datasheet</i>	Datum/Date: 19.06.2009 Ausgabe/Issue: 1 Datei/File: QSForm 116
--	---	--

Artikel	Druckinfusionsmanschette
<i>Article</i>	<i>Pressure Infusion Cuff</i>


REF	Mit Handgebläse, Dreiwegehahn und farbcodierter Druckanzeige <i>With hand inflator, 3-way stop cock and colour coded pressure indicator</i> 58-02-050 (500ml) 58-02-100 (1000ml) Mit Handgebläse, Dreiwegehahn und Präzisionsmanometer <i>With hand inflator, 3-way stop cock and precision manometer</i> 58-08-050 (500ml) 58-08-100 (1000ml)
------------	---

Beschreibung	Druckinfusionsmanschetten unterstützen die Verabreichung von Spül- und Infusionslösungen. Durch Belüften der Manschette wird Druck auf den Beutel ausgeübt.
<i>Description</i>	<i>Pressure Infusion Cuffs support the application of irrigation/infusion solutions. By inflating the cuff, pressure is exerted on the bag.</i>

Materialien (wenn vorhanden, Phthalat und dessen prozentualer Anteil			
<i>Materials</i> (if available, Phthalate and its percental rate)	1.	<i>Body:</i>	<i>Polyamide with TPU laminated</i>
	2.	<i>Gaze:</i>	<i>Polyester (Tetoron)</i>
	3.	<i>Hand Inflator:</i>	<i>Bulb: Soft PVC (Plastisol)</i> <i>Valve case: PVC</i> <i>Relief Valve: Rubber</i>
	4.	<i>Stop Cock:</i>	<i>Body: PVC</i> <i>Lever: HDPE</i>
	5.	<i>Pressure Indicator:</i>	<i>Body: ABS</i> <i>Inside magnetic spring</i>
	6.	<i>Tubing:</i>	<i>TPU</i>
	7.	<i>Connection:</i>	<i>Polyamide</i>
	8.	<i>Tubing:</i>	<i>Silicone</i>
	9.	<i>Hook:</i>	<i>TPU</i>
	10.	<i>Loop:</i>	<i>Nickel coated brass</i>
	11.	<i>Manometer:</i>	<i>ABS (case); PC (glas)</i>

REF 58-02-xxx

REF 58-08-xxx



Erstellt durch, Datum

Issued by, Date

H. Tucholl, 28.02.2019